

TP3: Spectrophotometric dosage

Objectives: .


Know the spectrophotometry technique.

- Know and use the relationship between absorbance and the concentration of an absorbing species in solution (Beer Lambert's law).
- Determine the concentration of a Methyl Orange solution from a calibration curve.

I-Operating principle of the spectrophotometer Spectrophotometry

is a quantitative analytical method which consists of measuring absorbance.

The spectrophotometer passes monochromatic radiation (of wavelength λ) through a cell of length L containing a colored solution. It then measures the absorbance A (a quantity linked to the quantity of light absorbed by the solution).




II- Beer Lambert's law ($A = \epsilon \cdot C \cdot L$)

The Beer-Lambert Law is a law that determines the relationship between absorbance and the concentration of an absorbing species in solution.

L = length of the tank (1 cm in general, with an accuracy of 1%) • C = concentration of the solution

- ϵ = Specific extinction coefficient which depends on the wavelength. ϵ also varies according to intermolecular forces and therefore the solvent used.



The graph representing absorbance as a function of concentration, called a **calibration line (or curve)**, allows the unknown concentration of a solution to be determined from the measurement of spectrophotometric absorbance.

The correlation coefficient, " R^2 ", and the "Y" axis at the origin of the regression line allow you to demonstrate the acceptability of the linearity of the data.

1-Preparation of Methyl Orange solutions

We have an aqueous solution of Methyl Orange S_0 with a mass concentration of $C_0 = 100 \text{ mg/L}$ ($3 \cdot 10^{-4} \text{ mol/L}$). We are trying to produce different solutions of Methyl Orange with concentrations of C_i to be determined by dilution (see table below).

Complete the four columns of the table.

No. solution n	Volume V_i of S_0 of water added (mL)	Volume Total (mL)	volume 50 mL	Concentration C_i of the diluted solution of Methyl orange (mg/L)	Absorbance A (unitless)
1				8 mg/L	
2			50 mL	10 mg/L	
3			50 mL	12 mg/L	
4			50 mL	15 mg/L	
5			50 mL	20 mg/L	
6			50 mL		


• Rinse a burette with Methyl Orange S_0 solution .

Then fill it with this same solution. • Prepare in

different 50 mL flasks the solutions S_1 to

S_6 :

- _ Pour the required volume of Methyl Orange solution (indicated in the table) then the corresponding volume of water,
- _ Complete the necessary volume with distilled water (up to 50mL), cap and shake.

**2- Measures**

The reading wavelength is constant: $\lambda = 465 \text{ nm}$.

For a reference (or "blank") solution set $A = 0$ (red "zero" push button).

_ Measure the absorbance of each of the solutions produced. _

Complete the last column of the table.

Data: Methyl orange: Its chemical formula is $\text{C}_{14}\text{H}_{14}\text{N}_3\text{O}_3\text{SNa}$
(Molar mass: 327.33 g/mol Solubility: 5.2 CH₃ g/L (at 20°C) and $\text{pK}_a = 3.4 \text{ CH}_3$)

