

Practical Work N°2

Simple Distillation of NaCl and KMnO₄ Solutions

1. Objectives

- Understand the principle of simple distillation.
- Separate a volatile solvent (water) from non-volatile solutes.
- Identify the ions present before and after distillation.
- Study the temperature evolution during phase change.

2. Theoretical Principle

Simple distillation allows separation of a volatile liquid from dissolved non-volatile substances.

When heating an aqueous solution:

- Water evaporates at 100 °C.
- Dissolved ions remain in the distillation flask.

Case 1: Sodium chloride



Case 2: Potassium permanganate



During boiling, temperature remains constant at the phase change (liquid → vapor).




Fig 1 : Simple Distillation Setup

3. Equipment and Reagents

Equipment:

- Round-bottom flask (250 mL)
- Distillation head
- Thermometer
- Liebig condenser
- Heating plate
- Beakers -Test tubes
- Bunsen burner (flame test)
- Retort stand and clamps
- Boiling chips

Reagents:

- Aqueous NaCl solution
- Aqueous KMnO₄ solution
- Silver nitrate solution (AgNO₃)
- Acidified Fe²⁺ solution
- Distilled water

4. Sample Preparation

- Prepare 100 mL of NaCl solution ($\approx 0.5 - 1$ M).
- Prepare 100 mL of KMnO_4 solution (dilute violet solution).

5. Experimental Protocol

Step 1: Ion detection before distillation

NaCl solution

- Add $\text{AgNO}_3 \rightarrow$ white precipitate (AgCl) \rightarrow presence of Cl^-
- Perform flame test ($\text{NaCl} +$ methanol) \rightarrow yellow flame \rightarrow presence of Na^+

KMnO_4 solution

- Observe purple color \rightarrow presence of MnO_4^-
- Add acidified $\text{Fe}^{2+} \rightarrow$ decolorization \rightarrow confirmation of MnO_4^-
- Flame test ($\text{KMnO}_4 +$ methanol) \rightarrow light purple flame \rightarrow presence of K^+

Step 2: Distillation Setup

- Introduce 100 mL of solution into the flask.
- Add boiling chips.
- Start water circulation in condenser.

Step 3: Heating and Data Collection

- Heat gradually.
- Record temperature every minute.
- Maintain a distillation rate of 1–2 drops/sec.
- Collect the distillate.

Step 4: Ion detection after distillation

Perform the same tests on the distillate:

- AgNO_3 test
- Fe^{2+} test
- Flame test
- Visual observation

6. Results

Expected observations:

- Distillate is colorless.

- No precipitate with AgNO_3 .
- No reaction with Fe^{2+} .
- No characteristic flame coloration.
- Temperature increases then stabilizes at $\sim 100^\circ\text{C}$.


Graph:

Temperature ($^\circ\text{C}$) vs Time (min)

→ Increasing curve followed by a plateau at boiling point.

7. Questions

1. Why is the distillate colorless?
2. complete the following reactions



3. Why do dissolved ions not evaporate?
4. What does the temperature plateau represent?
5. Why are boiling chips added?
6. What would happen if the solute were volatile?
7. What is the difference between simple and fractional distillation?
8. Can distillation remove all impurities from water? Explain.